Lactate VET test kit

For veterinary use only!

Veterinary test kit for quantitative in vitro determination of lactate in plasma with a solo analyser.

EUROLyser 

 Eurolyser Diagnostica GmbH Tel. +43 (0)662 / 43 21 00
Bayernstraße 11a Fax +43 (0)662 / 43 21 00-50
5020 Salzburg/Austria www.eurolyser.com

English

Order information

Order number: VT 0220
Order number: VT 0221
Order number: VT 2200

Indication

Lactate VET test kit
Lactate VET test kit
Lactate VET control kit

Kit size

16 tests/pack
6 tests/pack
1 x 5 ml (decision level)

Allow single test to warm up at room temperature (20 - 25 °C) by placing the test into the test kit rack for a minimum of 10 minutes before use. Put the test kit package back into refrigerator.

Summary

In the face of anaerobic metabolism lactate is synthesized from pyruvate. Muscle tissue is the predominant source of lactate production. Increased lactate concentration is a sensitive indication for a significant metabolic disturbance. Possible reasons are hypoxia and hypoperfusion. Diseases associated with increased lactate concentration are i.e. hypovolemic shock, local hypoperfusion (i.e. intestinal ischaemia in the face of ileus, aortic-thromboembolism), hypoxia, anaemia, or sepsis. In racehorses lactate is used to either check performance ability or i.e. exclude latent myopathy.

Method

Colorimetric test. Photometric endpoint measurement of the absorbance at 546 nm wavelength.

Measurement Range

0.1 – 15 mmol/l

When values exceed this range, samples should be diluted 1 + 1 with NaCl solution (9 g/l) and the results multiplied by 2.

Sample Material

Use 20 µl of plasma. Do not use serum!

For measurements right after taking the blood sample:
Use li-heparin samples and separate plasma from red blood cells immediately (within 5 minutes).

In case measurements cannot be done immediately:
Use Na-fluorid/EDTA samples, centrifuge and measure at the earliest.

Stability in plasma (after centrifugation and separation):
8 hours at 15 – 25 °C, 14 days at 2 – 8 °C.

Test Kit

ERS cuvette filled with mono reagent.

Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2 – 8 °C.
Opening the vacuum package may limit the reagent stability to 3 months (stored at 2 – 8 °C) from the date of opening. DO NOT FREEZE!

Warnings and Precautions

DO NOT INGEST! Avoid contact with skin and eyes. Observe all necessary precautions for the use of laboratory reagents.

Waste Management

Always refer to local legal requirements.

Reference Range

Canine: 0.5 – 3.0 mmol/l
Feline: < 1.0 mmol/l
Equine: 0.5 – 2.0 mmol/l

It is recommended that each laboratory establishes its own reference ranges.

Quality Control

For internal quality control the Eurolyser Lactate VET control kit is recommended. Order number: VT 2200

Precision

Reproducibility within-run:

Control high; N = 20; mean = 3.87 mmol/l; CV = 1.69%;
Control low; N = 20; mean = 1.49 mmol/l; CV = 2.10%;

Correlation

Equine sample correlation: N = 40

y (solo Lactate) = 1.0602 x (Cobas Lactate) - 0.2564;
 R^2 = 0.9577;

Interferences

No interference up to:

Ascorbic acid	30 mg/dl
Bilirubin	60 mg/dl
Triglyceride	2000 mg/dl
Haemoglobin	1000 mg/dl
Dopamine	10 mg/l
L-Dopamine	20 mg/l
Methyldopamine	10 mg/l
Glycolic acid	1200 mg/l

For further information on interfering substances refer to Young DS¹.

References

1. YOUNG, D.S. (2000) Effects of Drugs on Clinical Laboratory Tests. 5th ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press;

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Bestellnummer: VT 0220
Bestellnummer: VT 0221
Bestellnummer: VT 2200

Indication

Lactate VET test kit
Lactate VET test kit
Lactate VET control kit

Kit size

16 Tests/Packung
6 Tests/Packung
1 x 5 ml (decision level)

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English

ATTENTION!

Allow single test to warm up to room temperature (20 - 25 °C) for a minimum of 10 minutes before use!

- 1. Preparation of test system**
 - 1.1 Insert RFID card
 - 1.2 Place ERS cuvette into test kit rack
 - 1.3 Take off ERS cuvette cap
 - 1.4 Press „Measurement“ button and enter all required information using the touchscreen

-
- 2. Sample preparation**
 - 2.1 Aspirate 20 µl sample material from centrifuged sample tube

-
- 3. Test processing**
 - 3.1 Dispense 20 µl sample INTO THE LIQUID in the ERS cuvette
 - 3.2 Firmly apply cap onto ERS cuvette
 - 3.3 Place ERS cartridge into the solo analyser
 - 3.4 Start automatic sample processing by pressing the  start button on the solo analyser