

## Determine SMN1 and SMN2 copy number in a single reaction in less than four hours

The AmplideX PCR/CE SMN1/2 Kit provides a simple, sensitive, and accurate assessment of SMN1 and SMN2 copy number. It is optimized for use on widely established CE systems, requires minimal hands-on time, and involves just a single PCR reaction.






Figure 1 AmplideX PCR/CE SMN1/2 Kit\* Workflow







## AmplideX® PCR/CE SMN1/2 Kit: Target Product Specifications • SMN1 and SMN2 detection multiplexed in single PCR reaction • Similiar workflow to AmplideX® PCR/CE kits (FMR1, SMN1) **Assay Format** · Calibrator and control included • Copy number resolution up to at least 4 copies (0, 1, 2, 3, ≥4) Output/Performance • High degree of specificity for SMN1 and SMN2 detection Blood Sample Type Number of Reactions 50 Turnaround Time Less than 4 hours System Compatibility Applied Biosystems™ 3130, 3500, 3730 Series Genetic Analyzers March 2019 **Availability** A00001 Catalogue Number

\*Research Use Only. Not for use in diagnostic procedures.



Asuragen, Inc. 2150 Woodward Street, Suite 100 Austin TX 78744 asuragen.com 2000-071